



Contact Name: Leslie Mayeaux
Phone: 601-992-4121x206
Fax: 601-992-4645
E-mail: leslie@fingerpro.net

ACQUIRING FINGERPRO ID REMOTE WORKSTATION

FingerPro ID Remote Workstation for Mississippi healthcare facilities comes in two configurations. ***Livescan Device Configuration*** includes a certified livescan device for capturing fingerprints using inkless technology. ***Desktop Scanner Configuration*** utilizes a desktop scanner for capturing the fingerprint images from inked cards.

Acquiring FingerPro ID Remote Workstation with a livescan device

Step 1

Make sure that your facility complies with the requirements for the ***Livescan Device Configuration*** in the document “FingerPro ID Remote Workstation Installation Checklist”.

Step 2

AD&S will schedule your facility for the installation. The software will be remotely installed prior to your training. An AD&S representative will install the system and train your personnel.

Step 4

The facility will be responsible for the following: The cost of the live scan, service charges for installation and training time and FBI compression software (available by quote from AD&S).

Acquiring FingerPro ID Remote Workstation with a desktop scanner

Step 1

Make sure that your facility complies with the requirements for the ***Desktop Scanner Configuration*** in the document “FingerPro ID Remote Workstation Installation Checklist”.

Step 2

Purchase and send to AD&S, Inc. a certified desktop scanner for testing and configuration. Use the shipment address as it is specified below

Step 3

Your facility will be scheduled for the installation and your software will be remotely installed prior to your training.

Step 4

The facility will be responsible for the following: Service charges for installation and training time and FBI compression software (available by quote from AD&S)